

## Protocol 902\_16S PCR

### Material needed:

1. Bacterial genomic DNA (it is highly recommended to use the Qiagen Stool Kit according to our modified protocol and using the Retsch bead beater to prepare genomic DNA from complex mixtures of bacteria, such as intestinal content, or the quick and dirty method for bacterial isolates as described in protocol 901)
2. Primers:
  - fD1: AGA GTT TGA TCC TGG CTC AG
  - fD2: AGA GTT TGA TCA TGG CTC AG
  - rP1: ACG GTT ACC TTG TTA CGA CTT
  - working concentration is 10 $\mu$ M (stock concentration is 100 $\mu$ M)

### Procedure:

5 $\mu$ l	10X Taq polymerase buffer
1.5 $\mu$ l	MgCl <sub>2</sub>
1 $\mu$ l	10mM dNTP stock
1 $\mu$ l	primer 1 (of 10 $\mu$ M stock)
1 $\mu$ l	primer 2 (of 10 $\mu$ M stock)
1 $\mu$ l	primer 3 (of 10 $\mu$ M stock)
34.1 $\mu$ L	water
0.4 $\mu$ L	Qiagen Taq

Pipette 5 $\mu$ L of DNA sample into PCR tubes

### Program:

First step:	94°C	5 min
35 $\times$ repeat of:	94°C	1 min
	43°C	1 min
	72°C	2 min
last step:	72°C	7 min
	10°C	forever

### Reference:

1. Weisburg, W.G., S.M. Barns, D.A. Pelletier, and D.J. Lane. 1991. 16S ribosomal DNA amplification for phylogenetic study. *J Bacteriol* 173:697-703.